



Icahn School of Medicine at Mount Sinai LINCS Center for Drug Toxicity Signatures

Standard Operating Procedure: Quantitative RT-PCR PromoCell Cardiomyocytes

DToxS SOP Index:A-3.0

Last Revision: June 29, 2015

Written By: Evren U. Azeloglu

Approvals (Date): Joseph Goldfarb (6/30/15)
Marc Birtwistle (5/7/2015)
Eric Sobie (6/30/2015)
Ravi Iyengar (6/29/2015)

Quality Assurance/Control (QA/QC) steps are indicated with **green highlight**.

Metadata recording is highlighted with **yellow highlight and superscript indices**.

-
1. In order to assess the redifferentiation status of PromoCell cardiomyocytes (PCM), we compare the expression of the following cardiac differentiation markers using quantitative RT-PCR with mRNA isolated from undifferentiated PCM and differentiated PCM that were cultured at the same time. The cardiac markers are:
 - a. α -sarcomeric actinin (PCR primer from Qiagen, Cat: QT00082131)
 - b. Cardiac troponin T2 (PCR primer from Qiagen, Cat: QT00089782)
 - c. HERG channel (PCR primer from Qiagen, Cat: QT01003254)
 - d. Mef 2C (PCR primer from Qiagen, Cat: QT00053368)
 - e. Gata4 (PCR primer from Qiagen, Cat: QT00031997)
 - f. Nkx2.5 (PCR primer from Qiagen, Cat: QT00010619)
 - g. GAPDH, loading control (PCR primer from Qiagen, Cat: QT00079247)
 2. Measure the total RNA concentration for all the samples that are going to be assayed as collected by DToxS SOP A-1.0
 - a. We use a NanoDrop ND-1000 Spectrophotometer
 3. For each sample, mix the following contents from Superscript III First Strand Synthesis Kit (Invitrogen, Cat: 18080-051) in a thin walled 0.6 mL microfuge tube, ordered by descending volume:
 - a. 1 μ g of RNA (x μ L)
 - b. RNase-free water (8-x μ L)
 - c. Oligo(dT) (1 μ L)
 - d. dNTP mix (1 μ L)
 4. After each addition, mix by pipetting up and down, and quick spin on a tabletop mini centrifuge

5. Incubate mixture at 65°C for 5 minutes and move onto ice
6. Prepare the following cDNA synthesis mix by mixing the contents from the Superscript III Kit in descending volume order:
 - a. 10 µL mix for each sample made up of
 - i. 25 mM MgCl₂ (4 µL)
 - ii. 10x RT Buffer (2 µL)
 - iii. 0.1 M DTT (2 µL)
 - iv. RNase OUT (1 µL)
 - v. SuperScript III RT (1 µL)
7. Add the synthesis mix into the RNA solution, mix by pipetting up and down, and quick spin on a tabletop mini centrifuge
8. Incubate at 50°C for 50 minutes
9. Terminate reaction at 85°C for 5 minutes
10. Add 1 µL of RNase H and incubate at 37°C for 20 minutes to remove any leftover RNA
 - a. cDNA is stable at -20°C for up to six months
11. We use the Applied Biosystems 7500 Real-time PCR System with SYBR Green PCR Kit (Qiagen, Cat: 204143) and the primers listed in Step 1
 - a. Dilute each primer to 10 µM by adding RNase-free water according to the manufacturer's label
 - b. For each sample, add 6 µL of SYBR Green Master Mix, 1 µL primer, 1 µL cDNA, and 4 µL RNase-free water into a single well of a 96-well plate (Applied Biosystems, Cat: 4306737)
 - c. Using 7500 Software bundled with the instrument, quantify message levels using "standard experiment" setting via the crossing threshold metric.